


# NeXus scientific data format: Introduction and applications

**Pete R. Jemian, group leader**  
Beamline Controls & Data  
Acquisition Group  
Advanced Photon Source, USA



[www.nexusformat.org](http://www.nexusformat.org)  
[github.org/nexusformat](https://github.com/nexusformat)

**Könnecke *et al.*, *J Appl Cryst* 48.1 (2015) 301–305, 10.1107/S1600576714027575**

2017-08-25

This work was supported in part by the U.S. Department of Energy, Office of Science, under Contract No. DE-AC-02-06CH11357.

# A few things to know about NeXus

Uses HDF5 as file format: <https://www.hdfgroup.org>

- NeXus is a common data format for neutron, x-ray, and muon science.
- HDF5 is a hierarchical format:
  - components: groups, datasets, attributes, links
- NeXus
  - Defines the layout of information in an HDF5 file
  - Uses the HDF5 group structure to organize data
  - Provides a dictionary of common terms
  - NeXus *application definitions* can further define terms and layout specific to a scientific community for the exchange of data

## Basic NeXus structure

```
NeXus HDF5 data file
NXentry
  NXdata
    # datasets
  NXinstrument # optional
    # raw data structures
  NXsample # optional
    # sample metadata
  # other metadata
```

*J. Appl. Cryst.* (2015). **48**, 301-305  
<https://doi.org/10.1107/S1600576714027575>



# Introduction

- NeXus is an effort by an international group of scientists
- Started in 1995
- common data exchange format for neutron, X-ray, and muon experiments
- built on top of the scientific data format HDF5 ([hdfgroup.org](http://hdfgroup.org))
- Purposes
  - raw data : any relevant data associated with a scientific instrument or beamline
  - processed data : structures for exchange of data between applications

# *The Components of NeXus*

- Design Principles
  - group, field, attribute, link
- Data Storage Objects
  - base classes, application definitions, and contributed definitions
- Subroutines
  - NAPI, Python (and other) examples
- Scientific Community
  - you

# NeXus Design Principles

## Underlying HDF5 Data Model in each class


- Group
  - subdirectory
- Field
  - dataset
- Attribute
  - describes a group or field
- Link
  - a pointer

```
1  entry:NXentry
2      instrument:NXinstrument
3      detector:NXdetector
4      data:[]
5          @long_name = "strip detector 1-D array"
6      bins:[0, 1, 2, ... 1023]
7          @long_name = "bin index numbers"
8      sample:NXsample
9          name = "zeolite"
10     data:NXdata
11         @signal = "data"
12         @axes = ["bins", "bins"]
13         @bins_indices = [0, 1]
14         data --> /entry/instrument/detector/data
15         bins --> /entry/instrument/detector/bins
```

# NeXus Link

## like a pointer

- Avoids writing same data many times
- Use as needed to build NXdata (or other)
- Implemented as HDF5 Hard Link
  - @target attribute added to identify “original” location
- HDF5 External file links allowed
  - e.g., high-frame rate area detectors
  - Need to coordinate these files



# NeXus Data Storage Objects

Describes how the groups, fields, attributes, and links are to be used

## *base classes (54)*

- A set of components that are used to construct a data file
- Defines the set of terms that *might* be used in an instance of that class

## *application definitions (34)*

- Specifies a data structure for a given domain, such as small-angle scattering
- Defines the minimum set of terms that must be used

## *contributed definitions (10)*

- Propositions from the community for base classes or application definitions
- Other NXDL files for long-term archival by NeXus
- Considered as either in incubation or a special case not for general use



# NeXus base classes

## Different parts of a data file

- Set of terms that *might* be used in an instance of that class
- Can add additional terms
- Almost everything is optional
- Should cover most possible use cases
- Defines the spelling and meaning of such terms
- No single predefined place for *all* possible data
- Expect community contributions

## A few of the 54 base classes

### NXcrystal

A crystal monochromator or analyzer.

### NXdata

NXdata describes the plottable data and related dimension scales.

### NXdetector

A detector, detector bank, or multidetector.

### NXdetector\_group

Logical grouping of detector elements.

### NXdetector\_module

Geometry and logical description of a detector module.

### NXdisk\_chopper

A device blocking the beam in a temporal periodic pattern.

### NXentry

**(required)** NXentry describes the measurement.

### NXenvironment

Parameters for controlling external conditions

# application definitions

## Interface with analysis or ...

- Minimum set of terms that **must** be used in an instance of that class
- Can add additional terms
- Definition declares what is optional
- Should cover most possible use cases
- Defines the spelling and meaning of such terms
- No single predefined place for *all* possible data
- Expect community contributions

A few of the 34 application def.

### NXiqproc

Application definition for any  $I(Q)$  data.

### NXlauetof

This is the application definition for a TOF laue diffractometer

### NXmonopd

Monochromatic Neutron and X-Ray Powder diffractometer

### NXmx

functional application definition for macromolecular crystallography

### NXrefscan

This is an application definition for a monochromatic scanning reflectometer.

### NXreftof

This is an application definition for raw data from a TOF reflectometer.

### NXsas

raw, monochromatic 2-D SAS data with an area detector

# A very simple NXDL file


Each NeXus class is written in NeXus Definition Language (NXDL)

```
1  <?xml version="1.0" ?>
2  <definition
3    xmlns="http://definition.nexusformat.org/nxdl/3.1"
4    xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
5    xsi:schemaLocation="http://definition.nexusformat.org/nxdl/3.1 ../nxdl.xsd"
6    category="base"
7    name="verysimple"
8    version="1.0"
9    type="group" extends="NXobject">
10
11    <doc>
12      A very simple NeXus NXDL file
13    </doc>
14    <group type="NXentry">
15      <group type="NXdata">
16        <field name="counts" type="NX_INT" units="NX_UNITLESS">
17          <doc>counts recorded by detector</doc>
18        </field>
19        <field name="two_theta" type="NX_FLOAT" units="NX_ANGLE">
20          <doc>rotation angle of detector arm</doc>
21        </field>
22      </group>
23    </group>
24  </definition>
```

- NXDL file is XML
- Defines dictionary
  - Data type
  - Documentation
- Describes how data are related
- Marks data if item is required
- Each NXDL is versioned
- Governed by a Schema
- XML can be validated
- Continuous integration used to ensure integrity

# Data Units

- NeXus uses the UDUNITS standard (<http://www.unidata.ucar.edu/software/udunits>)
  - Very flexible
  - Allows mathematical expressions such as  $\text{g/cm}^3$
  - Supported by C library for programmatic conversion (also Python)
- Units appropriate to field: wavelength has NX\_WAVELENGTH type units
- Units are not presently validated
  - But, there are plans to check them using the C library code




# *Programs that read & write NeXus files*

see: <http://download.nexusformat.org/doc/html/utilities.html>

- Anything that can read or write HDF5 files
- Data Analysis applications:
  - At least 14 (including IDL, Matlab, PyMCA, and NeXpy)
- Tools from the HDF Group ([hdfgroup.org](http://hdfgroup.org)): HDFview, h5dump, ...
- Python package: h5py
- NeXus command-line utilities (from <https://github.com/nexusformat/code>)
- Validation of NeXus files: cnxvalidate & punx

# Example: 1-D data shown with NeXpy



# NeXus raw data file structure

## Standard dictionary to describe your instrument


- NXentry: data belonging to one run or scan
- NXinstrument: data that describes the instrument
- NXmonitor: incident intensity
- NXsample: all about the sample
- NXuser: all about the sample
- NXdata: (typically) links to plottable data
- Only NXentry and NXdata are required

```
HDF5 file
  NXentry
    NXinstrument
      NXsource
      NXcollimator
      NXattenuator
      NXdetector
    ...
    NXsample
    NXmonitor
    NXuser
    NXdata
  NXentry
```

# Example raw data tree

## LRMECS, IPNS, ANL: Chopper Spectrometer Pulsed Neutron Source

- Data collected in 2001
- For clarity here, no attributes shown






# Structure of a NeXus processed data file

root of HDF5 file		
<b>NXentry</b>		required
All data belonging to this processed data entry		
<b>NXprocess</b>	Data needed to describe this processing step	
	<b>input:NXparameter</b>	
	<b>output:NXparameter</b>	
	...	
<b>NXsample</b>		
<b>NXdata</b>		required
Links to data representing a default view		

# Example: EPICS Processed Data

## USAXS, APS, Python and h5py and Matplotlib



Wed Aug 23 19:59:59 2017


# Example: multi-dimensional data

## AF1410 steel SANS as multi-D

- The AF1410 SANS data might be re-written as a multi-dimensional data structure:
  - Sector angle (nuclear v. nuclear+magnetic scattering)
  - Aging time
  - $Q$
- The previous SANS data would need to be binned to the same  $Q$  steps
- Eliminates extra structure:
  - 1 NXentry group remains
  - 1 NXdata group remains
- Could include sector\_angle and/or aging\_time in the analysis model


```
cs_af1410_multi.h5 : NeXus data file
@creator = "hypothetical example"
@default = AF1410
AF1410:NXentry
  @default = AF1410
  AF1410:NXdata
    @signal = I
    @I_axes = Q
    @Q_indices = 2
    @aging_time_indices = 0
    @sector_angle_indices = 1
    I:NX_FLOAT64[10, 2, 77]
    Idev:NX_FLOAT64[10, 2, 77]
    Q:NX_FLOAT64[77]
    aging_time:NX_FLOAT64[10] = [0, .25,
      0.5, 1, 2, 5, 8, 10, 20 50, 100]
    @units = h
    sector_angle:NX_FLOAT64[2] = [0, 90]
    @units = degrees
    sasinstrument:NXinstrument
    sascollimation:NXcollimator
    sasdetection:NXdetector
    sassource:NXsource
    sassample:NXsample
```

# Example: 2-D data from SPEC hklmesh



## Example: EPICS Area Detector

Pilatus USAXS, APS, Dectris 1k, Python and h5py and Matplotlib



# Example: many samples

AF1410 steel

10 samples (aging series)

two (NIST) SANS each sample


- Steel sample, aging series
- Area detector image at each aging time
- Binned in horizontal & vertical sectors
- 1-D I(Q) have different lengths
- (lots of other metadata in example file)
- file: cs\_af1410.h5
- [https://github.com/canSAS-org/NXcanSAS\\_examples/tree/master/1d\\_standard](https://github.com/canSAS-org/NXcanSAS_examples/tree/master/1d_standard)

*Acta Metall* **41** (1993) 1869-1884

NeXus

```
cs_af1410.h5 : NeXus data file
@creator = xml2hdf5.py
@default = AF1410_10
AF1410_10:NXentry
  AF1410_a10:NXdata
    I:NX_FLOAT64[77]
    Idev:NX_FLOAT64[77]
    Q:NX_FLOAT64[77]
  AF1410_b10:NXdata
    I:NX_FLOAT64[76]
    Idev:NX_FLOAT64[76]
    Q:NX_FLOAT64[76]
  sasinstrument:NXinstrument
  sascollimation:NXcollimator
  sasdeteor:NXdetector
  sassource:NXsource
  sassample:NXsample
```

```
AF1410_1h:NXentry
  AF1410_a1h:NXdata
  AF1410_b1h:NXdata
AF1410_20:NXentry
AF1410_2h:NXentry
AF1410_50:NXentry
AF1410_5h:NXentry
AF1410_8h:NXentry
AF1410_cc:NXentry
AF1410_hf:NXentry
AF1410_qu:NXentry
```



# Example: multi-modal data

- NeXus structure is flexible
  - describe many sets of measurements
- Use data groups under a single entry group if from the same sample
- (Sample and wavelength data not shown for brevity)
- Or, use application definitions

root of HDF5 file		
<b>NXentry</b>		required
All data from multiple methods		
<b>NXinstrument</b>		
<b>NXsource</b>		
...		
<b>NXsample</b>		
<b>NXuser</b>		
<b>NXdata</b>		required
Links to data representing a default view		
...		
<b>NXsubentry</b>		
e.g. small-angle scattering application definition		
<b>NXsubentry</b>		
e.g. wide-angle scattering application definition		

## 2-D SAS/WAS images

```
sas-was-example.h5 : NeXus data file
@canSAS_class = "SASroot"
entry:NXentry
  @canSAS_class = "SASentry"
  sasdata:NXdata
    @canSAS_class = "SASdata"
    @I_axes = Q, Q
    @Q_indices = 0, 1
    @signal = I
    I:NX_FLOAT64[100,512]
    Qx:NX_FLOAT64[100,512]
    Qy:NX_FLOAT64[100,512]
    Qz:NX_FLOAT64[100,512]
  wasdata:NXdata
    @I_axes = Q, Q
    @Q_indices = 0, 1
    @signal = I
    I:NX_FLOAT64[256,256]
    Q_x:NX_FLOAT64[256,256]
    Q_y:NX_FLOAT64[256,256]
    Q_z:NX_FLOAT64[256,256]
```

# Examples: How to describe uncertainty

- These values are the *estimates* of uncertainty.
- By default, assumed as estimated standard deviation, but could specify as other
- Can be added to any field with an attribute:
  - Exact name of uncertainty dataset is flexible

## simple uncertainty

```
I:NX_FLOAT64[91]
  @units = 1/cm
  @uncertainties = Iesd
Iesd:NX_FLOAT64[91]
  @units = 1/cm
Q:NX_FLOAT64[91]
  @units = 1/nm
  @resolutions = Qdev
Qdev:NX_FLOAT64[91]
  @units = 1/nm
```

## uncertainty from multiple sources

```
I : float[m,n]
  @uncertainties=Idev
Idev : float[m,n]
  @components=I_uncertainties
I_uncertainties : (group)
  electronic : float[m,n]
    @basis="Johnson noise"
  counting_statistics : float[m,n]
    @basis="shot noise"
  secondary_standard : float[m,n]
    @basis="esd"
```



## Example: metadata

Rich metadata can be added, but is optional

```
sassample:NXsample
  @NX_class = NXsample
  @canSAS_name = sassample
  @canSAS_class = SASSample
  ID:NX_CHAR[34] = AF1410-10 (AF1410 steel aged 10 h)
  details:NX_CHAR[128] =
    transverse saturation magnetic field (1.6 T) applied in
    horizontal direction to clear magnetic domain scattering
```

```
sasinstrument:NXinstrument
  sasscollimation:NXcollimator
  sassdetector:NXdetector
  sassource:NXsource
    incident_wavelength:NX_FLOAT64 = 0.85
    @units = nm
    radiation:NX_CHAR[7] = neutron
    wavelength_spread:NX_FLOAT64 = 25.0
    @units = percent
```

## Example: metadata citation

- NeXus provides the *NXnote* to describe any additional freeform information not covered by the other base classes
- These tags and attributes were supplied *ad hoc*

```
sasnote:NXnote
citation:NXcollection
  @NX_class = NXcollection
  @canSAS_name = citation
  @tag = citation
  journal:NX_CHAR[11] = Acta Metall
    @tag = journal
  pages:NX_CHAR[9] = 1869-1884
    @tag = pages
  title:NX_CHAR[92] = Small-Angle Neutron Scattering
    Studies of Carbide Precipitation in
    Ultrahigh-Strength Steels
    @tag = title
  volume:NX_CHAR[2] = 41
    @tag = volume
  year:NX_CHAR[4] = 1993
    @tag = year
  authors:NXcollection
    @NX_class = NXcollection
    @canSAS_name = authors
    @tag = authors
    author_0:NX_CHAR[10] = A.J. Allen
      @tag = author
    author_1:NX_CHAR[11] = D. Gavillet
      @tag = author
    author_2:NX_CHAR[13] = J.R. Weertman
      @tag = author
```

# NeXus International Advisory Committee

<http://www.nexusformat.org/NIAC.html>

## Terms of Reference

The purpose of the NeXus International Advisory Committee is to supervise the development and maintenance of the NeXus common data format for neutron, x-ray, and muon science.

## How to contact the NIAC

Send an email to: [nexus-committee@nexusformat.org](mailto:nexus-committee@nexusformat.org)

## Membership

NIAC seeks a balanced representation of the international community. It will consist of at most one voting representative from each major neutron, synchrotron x-ray, and muon facility.

**The NIAC invites PSL-II to become a member.**

# Communications

## *Mailing lists*

See: <http://download.nexusformat.org/doc/html/maillinglist.html>

## *Bi-weekly teleconferences*

Discuss timely issues: <http://www.nexusformat.org/Teleconferences.html>

Next one: Tuesday, 12 September, 16:30 CEST (+0200)

## *NeXus Code Camp*

Approximately yearly, next one October 24-26, 2017, Diamond Light Source, UK

## *NIAC Meeting*

Every two years, next one 2018, location not decided yet

# Who uses NeXus?

Only a partial list from 2013



Add to this list:

- EPICS Area Detector HDF5 plugin
- Petra III at DESY, Hamburg, Germany

# Acknowledgement

- Mark Basham, Diamond Light Source, UK
- Herbert Bernstein, CIF (non-facility member)
- Aaron Brewster, Lawrence Berkeley Laboratory, USA
- Stuart Campbell, Brookhaven National Laboratory/NSLS-2, USA
- Bjorn Clausen, Los Alamos National Laboratory, USA
- Stephen Cottrell, Rutherford Appleton Laboratory, UK (Muon Representative)
- Ricardo Ferraz-Leal, SNS and HFIR at ORNL, USA
- Jens-Uwe Hoffmann, Helmholtz Zentrum Berlin, Germany
- Pete Jemian, Advanced Photon Source, USA (Documentation Release Manager)
- Mark Könnecke, Paul Scherrer Institut, Switzerland (Executive Secretary)
- David Männicke, Australian Nuclear Science and Technology Organisation, Australia
- Raymond Osborn, Argonne National Laboratory, USA (non-facility member)
- Tobias Richter, European Spallation Source, Sweden (Chair)
- Armando Sole, European Synchrotron Radiation Facility, France
- Jiro Suzuki, KEK, Japan
- Benjamin Watts, Swiss Light Source, Switzerland
- Eugen Wintersberger, DESY, Germany (Technical Manager)





Thank you for your attention.